

The advantage of
working
with Interchangeable Articulators



Accurate
t r a n s f e r



fast
workflow



Enable interchangeability of stone models between articulators of the same brand (when calibrated) providing convenience in handling and clinical communication x laboratory and digital equipment.

Many of the articulators on the market are not calibrated. In this way, stone models mounted on a given SAA cannot be transferred to another SAA with guaranteed accuracy, even though they are from the same brand. Bio-Art developed a standardization kit consisting of tools used to calibrate the articulators (models **A7Fix** and **A7Plus**). Using this resource is essential for the professional who has more than one articulator, as well as for use between dentists and dental technicians without the need to keep the models always in the same articulator. The calibration allows the interchanging of the models previously mounted on calibrated articulator plates.

The gauging device (calibration column) must be used periodically to check the accuracy between the articulators.



Extracted article from the book:

Logical - A clinical approach to occlusion